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Section 2 – Technical Manual (English)

This technical manual refers to the ReliefMed SmartPad medical device, developed by ReliefMed Group S.R.L., based in Lecce (Italy). The device is designed for the prevention and treatment of pressure ulcers in hospital and nursing home environments.

The device integrates:

- Automatic pressure detection system with sensors and ARM microcontroller.
- Integrated negative pressure therapy module (Class I during the first year of development).
- Protective polyethylene layer, high-resilience polyurethane foam, and therapeutic hydrogels (collagenase + hyaluronic acid, starting from the second year of design).
- Central digital control unit with user interface for healthcare professionals.

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0. Safety, Liability, and Intended Use

0.1 Purpose and Scope

This manual describes the maintenance, inspection, and management procedures for the ReliefMed SmartPad v1.0, intended for hospital, clinical, and home use under qualified medical supervision.

0.2 Intended Use

The ReliefMed SmartPad is designed for:

- **prevention and reduction of pressure ulcer risk in patients with limited mobility;**
- **therapeutic support in early-stage skin lesions through pressure modulation and microclimate control;**
- **integration with clinical data monitoring systems via ARM interface and biometric sensors.**

The device is not intended to replace medical care or pharmacological treatments.

0.3 Contraindications and Limitations

Use of the device is contraindicated in cases of:

- acute skin infections at the contact area;
- untreated burns or open wounds;
- presence of implantable devices interfering with negative pressure therapy.

Do not use the device under mechanical stress beyond design limits.

0.4 Manufacturer and User Responsibility

ReliefMed Group S.R.L. ensures compliance with safety and quality standards, provided that:

- the device is used according to this manual;
- maintenance is carried out only by authorized technical personnel;
- no structural or electronic modifications are made without approval.

Failure to follow these instructions or improper use voids the warranty and releases the manufacturer from all liability.

0.5 Disposal

At the end of its useful life (approx. 5 years or 2000 operating hours), the device must be disposed of as electronic and biomedical waste (EWC code 18 01 03*), following local regulations.

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English

1. Introduction and Purpose

This technical manual provides instructions for the proper routine and extraordinary maintenance of the **ReliefMed SmartPad**, designed for prevention and treatment of pressure ulcers using **Negative Pressure Wound Therapy (NPWT)** and **integrated sensor monitoring**.

2. Device Description and Components

The **ReliefMed SmartPad** consists of:

- **Polyurethane base layer** (load distribution and structural support)
- **High-density Polyethylene layer** (comfort and support)

- **Two medical-grade silicone layers** with micromodular cushions
- **Pressure, temperature, and humidity sensors**
- **ARM microcontroller** for real-time data processing
- **Negative pressure therapy module** with centrifugal pump
- **Hydrogel modules** containing collagenase and hyaluronic acid (*from year two of design*)

3. Routine Maintenance

Operations include:

- Periodic cleaning using a neutral, non-corrosive solution.
- Checking sensor and ARM module performance.
- Verifying electrical and data connections.
- Running diagnostics using ReliefMed Diagnostic Tool software.

4. Extraordinary Maintenance

Includes:

- Replacement or repair of damaged modules.
- Verification of polyethylene and polyurethane HR layers.
- Sensor recalibration through ARM interface.
- Firmware update of the electronic control unit.

5. Safety and Transport Procedures

- Keep the device in a horizontal position during transport.
- Avoid humidity and temperatures above 45°C.
- No specific thermal control required.
- Protect from shocks or excessive bending.

6. Technical Checklist and Service Log

Each maintenance operation must be recorded in the company logbook, including:

- Date of intervention
- Responsible technician
- Type of operation
- Verification results
- Replaced or updated parts

7. Illustrative Appendix

Includes technical CAD diagrams and conceptual schematics showing:

- Layer architecture
- ARM module and sensor layout
- NPWT airflow circuit

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© 2025 ReliefMed Group S.R.L. – Lecce, Italia